


MCP Supporter





MCP -> Phosphor Screen의 경우 Lead선이 뒤쪽 방향
(Guard Ring은 추가하는 것이 좋음)

Beam Medical Instrument 34036 대전 유성구 테크노10로 29 201호 Ph: +82-42-485-7112 Fax:+82-42-486-7112 EMAIL : kwonhae@bmii.co.kr			QUOTATION NO.	BM-190718-1
			DATE	2019. 07. 18
			사업자등록번호.	314-86-33855
				-
Origin	Korea	Term of validity	08/27/2019	
Supplier	BMI, Inc.	Shipping port		
Payment terms	Advance payment	Delivery terms		
Shipment date	2 months	Destination		
packing	standard packing	Offer No.		
REMARK.예원해(010-7674-7118)			Total amount (부가세포함) ₩2,540,000	
			대표이사 예원해 	

To : 고려대학교 담당자

Product No	Description	Qty(EA)	Unit Price	Amount(KRW)
1	빔프로파일 측정모듈 마운트 도면 참고	1	₩1,540,000	₩1,540,000
2	빔프로파일 측정모듈 및 마운트 설계비	1	₩1,000,000	₩1,000,000
		Amounts		₩2,540,000
		VAT	포함가	₩0

RFEA(retarding field energy analyzer)



1. MCP 앞 단에 4개의 Mesh를 둘 공간이 없음
2. RFEA를 위해서는 전류 측정이 필요한데 MCP 뒤에 Screen이기 때문에 전류 측정 불가능
-> 필요하다면 새로운 Line에서 사용하는 것을 고려해야 함