



20171124 Meeting

**Korea University
Lim. EunHun**





질소 Vent O-ring 교체
(AN계열 128)
여분 5개 Lab



질소 Vent O-ring 위치

-> 힘으로 넣을 시
O-ring 손상

KU Magnet



액체 헬륨 100L 주입



Table 3-2. Operating requirements for the magnet

LN2 Hold Time	Minimum of 7 days
LN2 Refill Volume	About 80 L
LN2 Yearly Consumption	About 4160 L
LHe Hold Time	Minimum of 100 days, expect 130
<u>LHe Refill Volume</u>	About 100 L
LHe Yearly Consumption	About 365 L

Thank You

